



The
Patent
Office

PCT/GB98 / 0 2 9 0 4
EAU

PRIORITY DOCUMENT

SUBMITTED OR TRANSMITTED IN
COMPLIANCE WITH RULE 17.1(a) OR (b)

09/509687

The Patent Office
Cardiff Road
Newport
Gwent
NP9 1RH

REC'D	1 3 NOV 1998
WIPO	PCT

I, the undersigned, being an officer duly authorised in accordance with Section 74(1) and (4) of the Deregulation and Contracting Out Act 1994, to sign and issue certificates on behalf of the Comptroller-General, hereby certify that annexed hereto is a true copy of the documents as originally filed in connection with the patent application identified therein.

In accordance with the Patents (Companies Re-registration) Rules 1982, if a company named in this certificate and any accompanying documents has re-registered under the Companies Act 1980 with the same name as that with which it was registered immediately before re-registration save for the substitution as, or inclusion as, the last part of the name of the words "public limited company" or their equivalents in Welsh, references to the name of the company in this certificate and any accompanying documents shall be treated as references to the name with which it is so re-registered.


In accordance with the rules, the words "public limited company" may be replaced by p.l.c., plc, P.L.C. or PLC.

Re-registration under the Companies Act does not constitute a new legal entity but merely subjects the company to certain additional company law rules.

Signed *AWB*

Dated *25/10/98*

tents Form 1/77

Patents Act 1977
(Rule 16)The
Patent
Office30SEP97 E306530-1 002638
P01/7700 25.00 - 9720693.2

Request for grant of a patent
 (Please read the notes on the back of this form. You can also get a statutory leaflet from the Patent Office to help you fill in this form)

The Patent Office

Cardiff Road
Newport
Gwent NP9 1RH

1. Your reference WJN/NV/WCM.63

2. Patent application number
(The Patent Office will fill in this part)

9720693.2 30 SEP 1997

3. Full name, address and postcode of the or of each applicant (underline all surnames)

University of Wales College of Medicine
Heath Park
Cardiff
CF4 4XN

Patents ADP number (if you know it)

-- 798546003

If the applicant is a corporate body, give the country/state of its incorporation

United Kingdom

4. Title of the invention

Assay

5. Name of your agent (if you have one)

Wynne-Jones, Laine and James

"Address for service" in the United Kingdom to which all correspondence should be sent (including the postcode)

Morgan Arcade Chambers
33 St Mary Street
Cardiff
CF1 2AB

Patents ADP number (if you know it)

1792002

6. If you are declaring priority from one or more earlier patent applications, give the country and the date of filing of the or of each of these earlier applications and (if you know it) the or each application number

Country

Priority application number
(if you know it)Date of filing
(day / month / year)

--

--

--

7. If this application is divided or otherwise derived from an earlier UK application, give the number and the filing date of the earlier application

Number of earlier application

Date of filing
(day / month / year)

--

--

8. Is a statement of inventorship and of right to grant of a patent required in support of this request? (Answer 'Yes' if:

YES

a) any applicant named in part 3 is not an inventor, or

b) there is an inventor who is not named as an applicant, or

c) any named applicant is a corporate body.

See note (d))

Patents Form 1/77

9. Enter the number of sheets for any of the following items you are filing with this form. Do not count copies of the same document

Continuation sheets of this form

Description 12

Claim(s) -

Abstract -

Drawing(s) -

10. If you are also filing any of the following, state how many against each item.

Priority documents -

Translations of priority documents -

Statement of inventorship and right to grant of a patent (Patents Form 7/77) -

Request for preliminary examination and search (Patents Form 9/77) -

Request for substantive examination (Patents Form 10/77) -

Any other documents (please specify) -

11. I/We request the grant of a patent on the basis of this application.

Signature  Date 30.9.97

Wynne-Jones, Laine and James

12. Name and daytime telephone number of person to contact in the United Kingdom

W.J. Newell - Tel 01222 229526

Warning

After an application for a patent has been filed, the Comptroller of the Patent Office will consider whether publication or communication of the invention should be prohibited or restricted under Section 22 of the Patents Act 1977. You will be informed if it is necessary to prohibit or restrict your invention in this way. Furthermore, if you live in the United Kingdom, Section 23 of the Patents Act 1977 stops you from applying for a patent abroad without first getting written permission from the Patent Office unless an application has been filed at least 6 weeks beforehand in the United Kingdom for a patent for the same invention and either no direction prohibiting publication or communication has been given, or any such direction has been revoked.

Notes

- If you need help to fill in this form or you have any questions, please contact the Patent Office on 0645 500505.
- Write your answers in capital letters using black ink or you may type them.
- If there is not enough space for all the relevant details on any part of this form, please continue on a separate sheet of paper and write "see continuation sheet" in the relevant part(s). Any continuation sheet should be attached to this form.
- If you have answered 'Yes' Patents Form 7/77 will need to be filed.
- Once you have filled in the form you must remember to sign and date it.
- For details of the fee and ways to pay please contact the Patent Office.

ASSAY

The present invention relates to an assay for measuring antibodies to the thyrotropin receptor and specifically to the use of a cell line transfected
5 with α -luciferase cDNA therein.

The thyrotropin receptor (hereinafter 'TSH-R') is known to regulate both the function and proliferation of the thyroid cell, and is stimulated by the hormone thyrotropin (TSH). The TSH-R is also a target for autoantibodies that either block the action of TSH (TBAb) (ie. act as hypostimulants or
10 inhibitors) or stimulate (hyperstimulate) the thyroid (TSAb).

Stimulation by TSAb is believed to be a mechanism operating in Graves' disease, whilst inhibition by TBAb is believed to be the case in idiopathic myxoedema. Patients with hyperthyroid Graves' disease produce antibodies which mimic the action of TSH, leading to chronic stimulation of
15 adenylate cyclase; and whereas the autoantibodies in some patients with idiopathic myxoedema are also able to bind to the TSH-R, nevertheless this does not result in an increase in intracellular cyclic adenosine 3',5'-monophosphate (cAMP).

Many assays have been developed to measure TSH-R autoantibodies.
20 The most widely used is a radioreceptor assay in which the binding of bovine 125 I-TSH to detergent-solubilised porcine TSH-R is inhibited by immunoglobulins or sera from patients suspected of having TSH-R autoantibodies. The main problems with this assay are that it uses a non-

human antigen; and that it measures binding and not biological activity, hence it is not able to distinguish between TSAb and TBAb.

It is clearly desirable to be able to distinguish between stimulation and inhibition, and therefore attempts have been made to develop a bioassay in which an effect of TSH (or its inhibition by TBAb) or TSAb is measured, such as the increase in cAMP. This may be done in thyroid slices or thyroid cells in culture and the greatest sensitivity, defined as the highest percentage of Graves' patients being positive, is achieved when the assay is performed in hypotonic medium. Both the radioreceptor and the bioassays are hampered by limited availability of biological material. To overcome this, a rat thyroid cell line (FRTL-5) has been developed, but here species differences may still be problematic. The recent cloning and sequencing of the TSH-R give unlimited access to recombinant human TSH-R.

As described by Libert *et al* in *Biochem. & Phys. Res. Comm.* 165 (3) 1250-1255 (1989) [all references herein are herein incorporated in their entirety, together with any cross-references therein], previous cloning of the dog thyrotropin receptor opened the way to molecular characterisation of the human TSH-R via isolation of human TSH-R cDNA clones; the analysis of the primary structure of the encoded polypeptide; and evidence that the recombinant molecule binds auto-antibodies found in patients with Graves' disease and idiopathic myxodema. The dog TSH-R cDNA (a 2.8 kb fragment) was used to hybridise a human thyroid DNA library. Sequencing of the resulting clones gave rise to a 2292 nucleotide residue open reading

frame encoding a 744 amino acid polypeptide having 90.3% convergence with the dog TSH-R₂. Transfection of the coding sequence in the pSVL vector of COS-7 cells allowed confirmation of the protein's ability to bind specifically TSH.

5 The cotransfection of Chinese hamster ovary (CHO) cells with a pSVL vector containing the coding region of human TSH-R led to the selection of cell lines particularly responsive to TSH or TSAB in terms of their cAMP accumulation (reported by Perret *et al* in *Biochem. & Biophys. Res. Comm.* 171 (3) 1044-1050 (1990)). Dose response curves of TSH-mediated cAMP
10 accumulation were reported for clones JP14, JP26 and JP28, and the number of receptors per cell were found to be highest in clones JP14 and JP09.

This work gave rise to the possibility of a bioassay in which cAMP production is measured in cells stably transfected with the human TSH-R in the presence of an autoantibody either alone (TSAB) or in the presence of
15 TSH (TBAB), (Ludgate *et al* in *Molec. & Cell. Endocrin.* 73 R13-R18). However, such an assay is not sufficiently robust for routine use since it takes several days to perform and requires tissue culture facilities.

Accordingly, the present invention provides an assay for TSH-R autoantibodies or TSH comprising:

- 20 (a) bringing into contact a test sample with cells from a clone expressing human TSH-R stably transfected with a reporter construct comprising cDNA of both a substrate capable of causing a measurable response when brought into contact with a corresponding protein and

a promoter containing cAMP response elements whereby increased cAMP levels increases expression of the enzyme;

(b) adding the corresponding substrate to cells thus contacted;

(c) measuring the response in the cells exposed to substrate; and

5 (d) comparing the response from test step (c) with the response from a standard or normal sample which has undergone steps (a) to (c).

Suitable reporter constructs are those in which the enzyme activity results in a colour change or emission of light. Examples of such enzymes
10 include chloramphenicol acetyl transferase (CAT), *Firefly* luciferase, *Renilla* luciferase, β -galactosidase, alkaline phosphatase, horse radish peroxidase or green fluorescent protein.

The cyclic AMP response element (CRE) of the reporter gene could be any promoter sequence or synthetic oligonucleotide which contains the CRE
15 consensus sequence, TGACGTCA, preferably as a number of tandem repeats. A suitable promoter is that for the glycoprotein hormone alpha subunit which contains tandem cAMP response elements, described by Kay *et al* in *Endocrinology* 134 (2) 568-573 (1994). Another example is a construct driving the CAT enzyme which has been described by Chatterjee *et*
20 *al* in *Molecular Endocrinology* 5(1) 100-109 (1991).

However, the assay according to the present invention preferably comprises, in step (a), the use of a luciferase cDNA driven by a promoter containing cAMP response elements; and, in step (b), the use of luciferin;

which means that the response measured in steps (c) and (d) is light output from the luciferinised cells. Preferably, the luciferase is Firefly luciferase although Renilla luciferase or the like would also be suitable.

Most preferably, the reporter construct is α -luciferase, being a
5 luciferase cDNA driven by the promoter for the glycoprotein hormone α
subunit, mentioned above. For example, the plasmid pA3luc may be
employed, having the glycoprotein hormone α subunit promoter introduced as
described by Maxwell *et al* in Biotechniques 7 276-30 (1989). The α -
luciferase is therefore 846 base pairs of 5' flanking sequence and 44 base
10 pairs of exon 1 of the glycoprotein hormone α subunit promoter in the plasmid
pA3luc. Alternatively, the CRE-containing sequence could be subcloned into
a commercially-available luciferase reporter system such as the pGEM-luc
vector from Promega. A further alternative is to use a plurality of plasmids
such as in the system available from Stratagene (CREB reporting system, no.
15 219010), which includes plasmids enabling a luciferase response to be
measured following an increase in cAMP.

The cells may be obtained as described in the reference mentioned
above by Perret *et al*. Alternatively, the human TSH-R could be subcloned
into any eukaryotic expression vector (of which pSVL is an example)
20 available from Stratagene, InVitrogen or the like. For selectivity, the more
recently-developed dual vectors which incorporate the antibiotic resistance
gene within the same plasmid, such as pcDNAIII (available from InVitrogen)
may be used. Otherwise, a separate plasmid for selection may be employed.

Preferably, the cells used in the assay (step (a)) are those identified as from clone JP09 in the above-mentioned reference which have been stably transfected within the order of 10^5 TSH-R per cell. More preferably, they are co-transfected with both α -luciferase cDNA and a puromycin resistance
5 encoding plasmid such as pSV₂Neo (available from Clontech) to allow selection of assay cells with puromycin. Surviving cells are then tested for luciferase activity in response to TSH.

Preferably, the cells would be lyophilised (freeze-dried) and provided in individual containers with one container being used per assay.

10 Another co-transfection may be carried out to provide the assay with a method of correcting for the number of cells seeded in a well during use, in the case where non-lyophilised cells are to be used. Since the Renilla luciferase construct is constitutive and has different substrate requirements from Firefly luciferase, it provides such a method. The same value would be
15 expected from every well whilst differences would reflect varying cell number. An appropriate plasmid for this transfection is the Renilla luciferase plasmid available from Promega, no. E2241. It contains the Herpes simplex virus thymidine kinase promoter upstream from Renilla luciferase, which is thus constitutively expressed.

20 Furthermore, when using intact rather than lyophilised cells, in order to prevent distortion of the assay by the presence of any TSH present in serum used in the cell culture medium, the serum should be charcoal-stripped at around 24 hours prior to assay.

The present invention therefore further provides a reporter construct comprising cDNA of both an enzyme capable of causing a measurable response when brought into contact with a corresponding protein and a promoter containing cAMP response elements whereby increased cAMP
5 levels increases expression of the enzyme, in particular wherein the enzyme is a luciferase.

Accordingly, the present invention further provides cells from a clone expressing human TSH-R stably transfected with a reporter construct comprising cDNA of both an enzyme capable of causing a measurable
10 response when brought into contact with a corresponding protein and a promoter containing CRE; the clone; and cDNA expressing the stably transfected human TSH-R.

Preferably, the assay according to this invention is carried out by means of a kit to enable fast and convenient results in a regular medicinal
15 biochemistry laboratory or hospital pathology or diagnostic laboratory. The present invention therefore further provides a kit for carrying out the bioluminescent assay of the present invention, which kit comprises:

- (a) cells from a clone expressing human TSH-R stably transfected with a reporter construct comprising cDNA of both an enzyme capable
20 of causing a measurable response when brought into contact with a corresponding protein and a promoter containing cAMP response elements whereby increased cAMP levels increases expression of the enzyme;

- (b) a standard sample for the assay;
- (c) medium for culturing and/or reconstituting the cells; and
- (d) instructions for carrying out the assay according to the present invention.

5 The corresponding protein and reagents relating thereto, and means for carrying out the response measurements may also be provided as part of the kit. For example, the kit may further comprise:

- (e) buffer for lysing the cells;
- (f) buffer for the reporter construct, preferably luciferase buffer;
- 10 and
- (g) corresponding protein, preferably luciferin, in buffer;
- and, optionally, a luminometer.

Alternatively, for example, in the case where luciferase/luciferin are employed in the assay, a separate, commercially-available kit may be employed such
15 as one of those available from the Promega Corporation. These commercially-available kits include no. E1483 wherein the luciferase is Firefly luciferase, and a dual-luciferase system no. E1910 which also employs Renilla luciferase.

The present invention will now be illustrated with reference to the
20 following non-limiting examples:

Example 1: Preparation of cells for use in luminescent assay for TSH and TSH-R antibodies - Use of Firefly luciferase

Chinese hamster ovary cells (CHO-K1), available from the ATCC number CCL61, were subjected to calcium phosphate co-transfection, using
5 standard protocols (Current Protocols in Molecular Biology, 1996, John Wiley & Sons Inc. section 9.1.4), with pSV₂Neo (available from Clontech) and a eukaryotic expression vector carrying the human TSH-R gene, pSVL-hTSHR (available from G. Vassart, IRIBIIN, Brussels, Belgium), as described by Perret *et al* (1990, *ibid*). The cells were selected with 400 µg/ml geneticin
10 (G418) and cloned by limiting dilution. Clone JP09 was isolated by this method (and is also available from Prof. G. Vassart), which expresses approximately 10^5 receptors per cell, as assessed in TSH binding experiments (described by Costagliola S, Swillens S, Niccoli P, Dumont J, Vassart G, Ludgate M in J Clin Endocrinol Metab 75 1540 *et seq.* (1992)).
15 Clone JP09 was co-transfected, again using the above-noted standard calcium phosphate method, with a eucaryotic expression vector carrying the puromycin resistance gene (available from, *inter alia*, InVitrogen, Stratagene etc.) and pA3Luc (available from V. Chatterjee, Univ. Cambridge, UK). Cells were selected in puromycin, 2.5 µg/ml and cloned by limiting dilution. Clones
20 were tested for light output (see Example 2) in response to bovine TSH (from Sigma, T8931). Clones giving a good TSH response were tested with a panel of normal human sera, the selection criterion in this case being a low light output of ≤ 1.5 relative light units.

Example 2: Preparation of cells for use in luminescent assay for TSH and TSH-R antibodies - Use of Firefly and Renilla luciferases

The method of Example 1 is followed, but after the second co-transfection involving the puromycin resistance gene, another co-transfection is carried out repeating the standard calcium phosphate method, but instead with a eukaryotic expression vector carrying the hygromycin resistance gene (available from Stratagene) and the R. luciferase plasmid E2241 (available from Promega).

Example 3: Assay for Measurement of (i) TSH bioactivity and (ii) Thyroid Stimulating (TSAB) or (iii) Blocking (TBAB) Antibodies

Culture medium for the cells is RPMI: 10% foetal calf serum (FCS); 1% glutamine; 1% pyruvate; 2% penicillin/streptomycin (and 2.5µg/ml puromycin, when amplifying cells for assay). 96 well plates are seeded with 5×10^4 cells prepared according to Example 1 and cultured overnight (approximately 16 hours) at 37°C in a water-saturated incubator. They are then cultured for a second night in medium containing 10% charcoal stripped calf serum (available from Gibco). Basal light output is measured, in triplicate wells, in the presence of 100µl of the medium containing 10% charcoal stripped serum. Standards (see below) and test samples are in the form of 10% patient serum to a final volume of 100µl of the medium containing 10% charcoal stripped calf serum. Basal, standard and test incubations are for 2-3 hours as above. The final result is expressed as relative light units (RLU) obtained by the ratio of standard:basal or test:basal.

(ii) TSH bioactivity: this test is for patients having high circulating TSH levels when measured by radio-immunoassay but low circulating free T4, indicative of hypothyroidism because of defective TSH. Standards are well-characterised euthyroid and hypothyroid serum samples of increasing biological activity.

(iii) TSAB measurements are made in thyrotoxic patients. Standards are pooled normal human serum (having a light output of ≤ 1.5 RLU), bTSH and well-characterised TSAB containing sera of increasing activity.

(iii) TBAB measurements are made in the same patients as in (i), but the assay wells also contain 1mU/ml of bTSH. Standards are pooled normal sera + 1mU/ml bTSH and well-characterised TBAB containing sera of increasing activity.

Measurement of light output

Following the incubation period, supernatants are removed from the wells and the cells are washed twice in phosphate buffered saline. Light output is measured using a commercially available kit from the Promega Corporation, no. E1483, according to the manufacturer's instructions. This involves treating the cells in lysis buffer (Promega no. E1513), adding the Firefly luciferase reagent and measuring the light output in a luminometer.

If the cells used for the assay also express Renilla luciferase constitutively (as prepared in Example 2), to give a method of standardisation of cell number/well, the Promega Dual-Luciferase system is used, no. E1910.

In this case, following measurement of the Firefly luciferase as above, a reagent to quench the luminescent signal is added followed by the Renilla luciferase reagent and a second reading taken.